



BIOTECK GRAFTS & MEDICAL DEVICES FOR REGENERATIVE MEDICINE

FOR INTERNAL USE ONLY


Fabrizio Raimondi
Business Development Manager

BIOTECK S.p.A.


BIOTECK is an Italian Company, manufacturing horse-derived substitutes for bone, cartilage and soft tissues reconstruction in Orthopaedics, Spine, General and Oral-Maxillofacial Surgery. Founded in 1995, it has developed Zymo-Teck[®], a proprietary enzymatic de-antigenation process, which guarantees grafts with optimal biological and biomechanical properties, in a wide range of formats (chips, blocks, paste, putty, crunch, membranes).

OVER 20 YEARS OF TRUSTABILITY

Founded in 1995, it has developed Zymo-Teck[®], a proprietary enzymatic de-antigenation process, which guarantees grafts with optimal biological and biomechanical properties, in a wide range of formats (chips, blocks, paste, putty, crunch, membranes).



Commercial headquarter in Vicenza



Production facility,
R&D and Customer Center

PRESENCE IN FOREIGN MARKETS



ARGENTINA
CHILE
ECUADOR
COLOMBIA
BRASILE
PANAMA
COSTARICA
REPUBBLICA DOMINICANA
SAN SALVADOR
GUATEMALA
HONDURAS
MEXICO
SPAGNA
PORTOGALLO
FRANCIA
OLANDA
GRAN BRETAGNA
IRLANDA
DANIMARCA
GERMANIA
SVIZZERA
AUSTRIA
SLOVENIA

SERBIA
GRECIA
ALBANIA
KOSOVO
BULGARIA
ROMANIA
POLONIA
SLOVACCHIA
REPUBBLICA CECA
BIELORUSSIA
UNGHERIA
UCRAINA
RUSSIA
ESTONIA
LITUANIA
LETTONIA
ARMENIA
AZERBAIJAN
IRAN
TURCHIA
UAE
ARABIA SAUDITA
EGITTO

QATAR
BAHRAIN
SYRIA
KUWAIT
CIPRO
GEORGIA
LIBANO
GIORDANIA
PALESTINA
IRAQ
KIRGHIZISTAN
HONG KONG
INDIA
THAILANDIA
MALAYSIA
SINGAPORE
GEORGIA
AFGHANISTAN
PAKISTAN
LIBANO
GIORDANIA
GEORGIA
AFGHANISTAN

PAKISTAN
LIBANO
GIORDANIA
PALESTINA
IRAQ
KIRGHIZISTAN
HONG KONG
INDIA
THAILANDIA
MALAYSIA
SINGAPORE
SINGAPORE
INDONESIA
FILIPPINE
TAIWAN
KOREA DEL SUD
REUNION
POLINESIA
FRANCESE
SUDAFRICA
ALGERIA
SUDAN

BIOTECK bone grafts unique features

Bone substitutes, manufactured by Bioteck SpA are natural grafts derived from heterologous equine bone tissue subjected to a proprietary patented cleaning process named Zymo-Teck®.

This tissue processing allows achieving fully biocompatible implants while preserving all their biological and biomechanical features for an optimal healing process after their implantation.

Zymo -Teck® is a multi-step process that guarantees the removal of all potentially antigenic and/or immunogenic components from the treated bone, by specific enzyme mixtures effective against specific molecular targets such as lipoproteins, glycoproteins and collagen telopeptides.

Several baths in mixtures of glycolytic and lipolytic enzymes take place with pressurized washing cycles with osmotic water at low controlled temperature. All products are then treated through an oxidative phase using hydrogen peroxide for the removal of cellular debris and any other contaminants. Finally, grafts are freeze-dried, packed in double blister and sent to β -ray terminal sterilization at 25kGy.

This type of irradiation are more respectful of biological and bio-mechanical behaviors of grafts than the γ rays and is allowed only due to the very low bio- burden (level of contamination) of Bioteck substitutes at the end of the tissue processing.

The Zymo -Teck® process does not use high temperature and chemical solvents. This ensures the complete preservation without any even minimal alteration of bone collagen and mineral component, thus making them pure biological matrix acting as natural scaffold (optimal osteo-conduction capability) to newly-bone formation. The presence of bone collagen in its intact structure allows Bioteck grafts to properly interact with cellular elements involved in the healing process, favoring the formation of new living tissue in physiological time and manner.

Bone collagen in its native structure exerts all the effects ascribed to it and, in particular , interacts with the sub-units of the beta-1 integrin of osteoblast membrane enhancing cell adhesion to scaffold¹ ; It acts as a co- activator of morphogenic proteins increasing the stimulating action of local growth factors² ; it interacts with mesenchymal cells from bone marrow by inducing their adhesion, proliferation and osteoblastic differentiation³ and finally, when implanted in a bone defect , exerts a pro - regenerative direct promoting the neo- osteogenesis⁴ .

It is also important to consider that the preservation of endogenous collagen structure , not only gives important biological properties, but it also helps to make them more resilient and resistant with excellent strength and load-bearing properties, particularly required in many orthopedic indications.


Bioteck grafts must be rehydrated in saline solution according to the instructions, before implantation. It is highly recommended to use when possible AWAYR perfusion device for a complete saturation of Bioteck porous bone grafts with saline solution or patient own biological fluids, like bone marrow, peripheral blood, PRP, etc. in order to remove all the air inside grafts while forcing cells and growth factors inside its structure thus enhancing and speeding up bone integration and ingrowth.

Complete antigens removal and total compatibility profile of Bioteck bone substitutes make them the safest and most reliable options for bone defect repair or regeneration. Hundreds of thousands of implanted grafts worldwide, with no report of adverse events together to many clinical experiences and scientific studies demonstrate that Bioteck grafts are a suitable material in many surgical applications even in presence of severe bone defects.

References

- 1) Type I collagen in xenogenic bone material regulates attachment and spreading of osteoblasts over the Beta 1 integrin subunit
Baslè, Lesourd, Grizon, Pascaretti, Chappard
Orthopade 1988 Feb. 27 (2) 136-42
- 2) Dissociative extraction and reconstruction of extra-cellular matrix components involved in local bone differentiation
Sampath, Reddi
PNAS 1981 Dec; 78 (12) 7599-803
- 3) Effect of type 1 collagen on the adhesion, proliferation and osteoblastic gene expression of bone marrow derived mesenchymal stem cells
Liu, Hu, Zhao, Wu, Xiong, Lu
Chin J Traumatol. 2004 Dec. 7 (6) 358-62
- 4) Evaluation of the effect of heterologous type 1 collagen on healing of bone defects
Gungormus, Kaya
J Oral Maxillofacial Surg. May. 60 (5); 541 -S

COMPOSITION OF BONE TISSUE




Antigens: Proteins, Polysaccharides, Glycoproteins, Lipoproteins and LPS


Antigenicity: interaction with an antibody

Immunogenicity: induction of an immune response

THE ZYMO-TECK[®] MOLECULAR TARGETS



WHY EQUINE TISSUE



It is considered one of the safest animals.

No prions pathologies, non transmissible diseases.

Not included in 2003/32/EEC Directive

Bred in fence – good bone morphology

SAFETY VS QUALITY

In the production of biological grafts SAFETY and QUALITY are strictly connected and directly depending on the applied process.

Most of the time in the absence of a specific biotechnology, in order to guaranty the safety many biological features are sacrificed.



THE ZYMO-TECK® PROCESS STEPS



5 - LYOPHILIZATION AND PACKAGING



Freeze dry Lyophilization


*ISO 7 Clean Room
and ISO 5 Laminar Flow hood*



Double blister packaging


Once packed in their final box, items are sent to terminal sterilization at 25kGy β rays.

QUALITY CONTROLS



ENVIRONMENTAL CONTROLS

Bioteck continuously checks the *Microbiological Quality*



Bioteck Microbiology Laboratory

STANDARD GRAFTS

Bioteck offers a wide range of grafts with over 300 different available codes.

Standard grafts belong to the OSTEOPANT® line, made of osteoconductive cancellous, cortical and cancellous-cortical grafts.

Load bearing capability for cancellous blocks have been tested higher than 360 kg/cm²

(Mechanical tests performed by University of Padova)



Cancellous Chips



Cancellous Blocks



Cancellous-cortical Strips



Cortical Rod



Cancellous Dowels



Hemi-femoral Head



HTO Wedge for plate fixation



HTO Wedge

BIOTECK
The science of bone tissue

safety/biocompatibility

osteoconduction

complete remodeling

naturalscaffold

osteopant
osteoplastflex
biocollagen

cancellous bone chips		
06-01-05	Cancellous Chips	(4 - 6 mm) 5 cc
06-01-10	Cancellous Chips	(4 - 6 mm) 10 cc
06-01-20	Cancellous Chips	(4 - 6 mm) 20 cc
06-01-30	Cancellous Chips	(4 - 6 mm) 30 cc
06-01-40	Cancellous Chips	(4 - 6 mm) 40 cc
06-01-60	Cancellous Chips	(4 - 6 mm) 60 cc

cancellous blocks		
09-01	Cancellous Block	30 x 20 x 30 mm
09-01A	Cancellous Block	10 x 15 x 30 mm
09-01B	Cancellous Block	10 x 10 x 20 mm
09-01C	Cancellous Block	10 x 10 x 30 mm / 2 pc
09-02	Cancellous Block	50 x 40 x 5 mm
09-02B	Cancellous Block	40 x 30 x 30 mm
09-03	Cancellous Block	50 x 40 x 30 mm

bio-gen putty		
84P-02	810-12K Putty	2 cc
84P-05	810-12K Putty	5 cc
84P-30	810-12K Putty	5 cc / 2 pc

cancellous dihedron		
09-07A	Cancellous Dihedron	50 x 30 x 30 mm

cancellous wedges		
09P-05	Cancellous Wedge	40 x 30 x 30 mm
09P-05B	Cancellous Wedge	40 x 30 x 30 mm
09P-06	Cancellous Wedge	50 x 40 x 30 mm
09P-06B	Cancellous Wedge	50 x 40 x 30 mm
09P-07	Cancellous Wedge	50 x 20 x 30 mm

cancellous wedges for Plating Fixation		
09P-01SP	Conc. Wedge for Plating Fixation	50 x 40 x 2,5 mm
09P-01SP	Conc. Wedge for Plating Fixation	50 x 40 x 3,5 mm
09P-012SP	Conc. Wedge for Plating Fixation	50 x 40 x 12,5 mm

flex acetabular-mat		
09F-070	Flexibla Acetabular-Mat	ø 70 x 5-7 mm

flex cancellous sheets		
07C-02	Flexibla Cancellous Sheet	40 x 20 x 3 mm
07C-04	Flexibla Cancellous Sheet	30 x 20 x 3 mm
07C-04	Flexibla Cancellous Sheet	50 x 25 x 3 mm
07C-05	Flexibla Cancellous Sheet	50 x 40 x 3 mm

flex cortical sheets		
07C-04	Flexibla Cortical Sheet	40 x 40 x 1,2-3 mm
07C-05	Flexibla Cortical Sheet	50 x 25 x 1,2-3 mm
07C-07	Flexibla Cortical Sheet	50 x 50 x 1,2-3 mm
07C-08	Flexibla Cortical Sheet	70 x 70 x 1,2-3 mm
07C-08	Flexibla Cortical Sheet	40 x 40 x 0,7-1 mm


biocollagen fleeces		
80G-005	Biocollagen Fleece	75 x 50 x 8 mm
80G-008	Biocollagen Fleece	50 x 80 x 8 mm
80G-100F	Biocollagen Fleece	300 x 80 x 8 mm

biocollagen membrane		
80G-07	Collagen Membrane	70 x 50 x 0,2 mm

hemi-femoral head		
09F-04	Hemi-Femoral Head	ø 60 mm

cortical plates		
09P-08	Cortical Plate	60 x 20 x 6 mm
09P-09	Cortical Plate	300 x 20 x 6 mm
09P-10	Cortical Plate	120 x 10 x 6 mm
09P-12	Cortical Plate	100/200 x 20 x 6 mm

FLEXIBLE GRAFTS



The OSTEOPANT Flex grafts are partially demineralized through a process of hydrolysis of calcium salts in acid environment that makes bone collagen more exposed compared to a standard graft.

For this reason, OSTEOPANT® Flex grafts promote cell adhesion and reduces the time of incorporation and remodeling.



Osteopant Cortical Flex Sheet



Osteopant Cancellous Flex Sheet



Osteopant Acetabular Mat



Osteopant Cancellous Flex Disc

XENO-DBM BONE PASTES

Cortical bone subjected to a treatment with selective enzymes and completely demineralized. Activagen contains the molecular signals typical of the extra-cellular matrix. These granules are the basis for our osteoproduktive bone paste.



Osteplant Activagen Injectable Paste



Osteplant Activagen Mouldable Paste.

<p>BIOTECK The science of bone tissue</p> <p>safety/biocompatibility</p> <p>osteopromotion</p> <p>boneregeneration</p> <p>bonepaste biocollagen gel & crunch osteoplant activagen osteoplant angiostad</p>	<p>biocollagen gel</p> <p>It is an injectable gelatine of type I and III collagen extracted from equine Achilles tendon, combined with cancellous bone powder (0.4 mm) to increase its consistency. It may be used for concrete grafts, thus promoting their stabilization and having an haemostatic action so it may be used as a carrier for granular grafts, platelet growth factor and autologous bone marrow concentrate.</p> <p>BSG-GBL2 Biocollagen Gel syringe - 7cc BSG-GBL5 Biocollagen Gel syringe - 5cc BSG-GBL10 Biocollagen Gel syringe - 10cc</p>		<p>osteoplant activagen injectable paste</p> <p>It is an injectable bone paste based on demineralized bone matrix of equine cortical bone in a collagenous carrier. It contains the DBM-specific molecular signals that promote the cascade process of bone regeneration. It may be used alone in traumatic injury in comminuted fracture and slight cavitary bone defects, or combined with other grafts, platelet growth factors and autologous bone marrow concentrate.</p> <p>OSG-AC15 Osteoplant Activagen Injectable Paste - 5cc OSG-AC10 Osteoplant Activagen Injectable Paste - 10cc</p>	
	<p>biocollagen crunch</p> <p>It is an osteoconductive bone paste based on type I and III collagen extracted from equine Achilles tendon, combined with bone powder and cancellous chips (0.4-2 mm), in ready-to-use syringes. It is recommended in case of filling bone defects, alone or with autologous bone grafts, platelet growth factors and autologous bone marrow concentrate. It is very mouldable and has good properties of site adhesion. Furthermore, the great quantity of collagen fosters the blood clot formation.</p> <p>BSG-CR15 Biocollagen Crunch syringe - 5cc BSG-CR10 Biocollagen Crunch syringe - 10cc</p>		<p>osteoplant activagen mouldable paste</p> <p>It is a mouldable bone paste based on demineralized bone matrix of equine cortical bone and cancellous bone chips (0.4-2 mm) in a collagenous carrier, packed in cut-off syringes. It contains the DBM-specific molecular signals that promote the cascade process of bone regeneration. It may be used alone or combined with platelet growth factors and autologous bone marrow concentrate.</p> <p>OSG-AM10 Osteoplant Activagen Mouldable Paste syringe - 0.5cc OSG-AM1 Osteoplant Activagen Mouldable Paste syringe - 1cc OSG-AM2 Osteoplant Activagen Mouldable Paste syringe - 2cc OSG-AM5 Osteoplant Activagen Mouldable Paste syringe - 5cc OSG-AM10 Osteoplant Activagen Mouldable Paste syringe - 10cc</p>	
	<p>osteoplant angiostad</p> <p>Angiostad is an injectable gel specifically aimed at promoting neoangiogenesis. It is composed of a demineralized matrix containing signals which support the formation of new capillaries thus fostering graft vascularization. It is recommended in those cases in which regeneration may be difficult (i.e. when the site between the volume to be regenerated and the vital bone surface is adverse).</p> <p>OSG-AB2 Osteoplant Angiostad syringe - 2cc</p>		<p>BIOTECK BIOTECK</p> <p>BIOTECK bone pastes are ready for use and can be prepared at room temperature. They are given during the bone production phase. Due to their biological characteristics and their great safety, they are recommended. In difficult situations, the use of some bone matrix, collagenase, bone granules, autologous serum, platelet growth factor and autologous bone marrow concentrate increases the great quantity of type I collagen.</p> <p>BIOTECK bone pastes on a class II Medical Device.</p> <p>CE 0173</p>	

Animal DBM works in human

Bone morphogenetic proteins and growth factors have been highly conserved during the evolution. This translates into a very high level of correspondence between human BMPs and those of other mammals. Since many years there are in literature positive results on the effectiveness of DBM from different species.

BMPs are not species-specific. High level of identity of amino-acid sequence of BMP in different animal species (Wozney JM, Cellular and Molecular Biology of Bone 1993)

Evolution of the transforming growth factor-beta superfamily. - Burt DW, Law AS. Prog Growth Factor Res. 1994;5(1):99-118.

Homology of bone-inductive proteins from human, monkey, bovine, and rat extracellular matrix. Sampath TK, Reddi AH. Proc Natl Acad Sci U S A. 1983 Nov;80(21):6591-5

Implants of heterologous (bovine) demineralized bone matrix for induction of posterior spinal fusion in rats. Spine. Guizzardi S, Di Silvestre M, Scandroglio R, Ruggeri A, Savini R. 1992 Jun;17(6):701-7

Healing response to various forms of human demineralized bone matrix in athymic rat cranial defects. Chesmel KD, Branger J, Wertheim H, Scarborough N. J Oral Maxillofac Surg. 1998 Jul;56(7):857-63; discussion 864-5

ACTIVABONE® - NEXT GENERATION OF NATURAL BONE PASTE



ACTIVABONE® is the innovative line of heterologous bone pastes with high biological activity and superior handling properties, based on the new technological platform of **Bioteck®** proprietary hydrogel with bio-modulated viscosity (Exur®).

This new carrier gives bone pastes better chemical-physical and rheological performances, allowing for suitable in situ stability with proper resistance to bleeding and irrigation.

The development of hydrogel-based injectable, mouldable and pre-formed **ACTIVABONE®** bone pastes has been foreseen as practical way to also enhance adaptability to defects of irregular geometry maintaining full contact with the live bone surrounding the defect area. The novel formulations guarantee longer maintenance of osteoconductive granules and osteoinductive DBM powders into defect and favour better attachment of osteogenic cells and growth factors to graft. These excellent biological and physical properties therefore translate into an optimal bone defect regeneration, by promoting and speeding-up early phases of tissue healing process with an improvement of the graft integration and remodeling.

This line of new bone pastes will be available on the market by the end of 2016.

EXUR® Hydrogel Proprietary Technology

For the past three years, Bioteck® has been focusing its effort and attention to medical grade polymers commonly used in pharmaceutical applications, such as poly(ethylene glycol) - PEG, poly(ethylene oxide) – PEO and hydroxypropyl methyl cellulose – HPMC, and developed an innovative proprietary technology to control their sinergic polymerization.

The ability to regulate polymerization and cross-linking density of polymers avoids granules dispersion and loss during surgery, for assuring complete filling and direct contact with the tissue surrounding the defect and successful bone repair. HPMC, PEG and the chemically similar PEO hydrogels undergo a polymerization reaction by physical sterilization, which is modulated by introducing very limited amount of Vitamin C. This anti-oxidant molecule is able to limit intra- and inter-molecular rearrangement of PEG and HPMC polymeric chains, then specifically tailoring visco-elasticity of Exur®, the **Bioteck®** proprietary hydrogel carrier. The possibility to select specific polymers molecular weight and concentration, as well as finely tuned Vitamin C amounts, therefore allows **Bioteck®** to design unique bone fillers with different physical and handling properties, either injectable, mouldable or shapeable pre-formed bone substitutes.

COLLAGEN FLEECES & MEMBRANES

Biocollagen Membrane & fleeces are medical devices mainly composed of type I collagen from equine achilles tendon. Fleeces are used as graft extensor together to bone chips or for haemostatic use. Collagen membrane main indications are the protection of the grafted site or periosteum substitute.



Biocollagen Fleece



Biocollagen membrane



BIOCOLLAGEN MeRG



It is a particular membrane for cartilage treatment derived from Biocollagen. It is a dual side collagen membrane indicated in the treatment of chondral lesions with AMIC technique

(Autologous Membrane Induced Chondrogenesis)

The rough and the smooth side at the SEM



MERG: DESCRIPTION

Re-absorbable membrane for guided tissue regeneration, (Biocollagen® MeRG).

Composition:

Mainly composed of Type I non-allergenic lyophilised collagen from equine Achilles' tendon (Type I >95%. See the following slides for specific instrumental analysis on contents.

Properties:

Microfibrillar collagen membrane. The peculiar arrangement of its fibers, given by the exclusive manufacturing process, gives the membrane resistance to torsion, pulling and tearing. The membrane features a smooth and a rough/fibrillar side.

The rough side improves adhesion, and allows for fixation with fibrin glue. Adhesion increases if blood is present. The collagen the membrane is made of is achieved from equine Achilles' tendons (one of the principal sources of type I collagen). Equines can not transmit encephalopathies.

Collagen shows the chemical and structural features of a glycoprotein, capable of interacting with the fibroblasts' and platelets' receptors. It activates factors XII and VIII and constitutes the structural basis of connective tissues. Its interactions with platelets is fundamental for coagulation, since its bonding with the platelets integrins leads to platelet degranulation and to release the factors that activate coagulation.

This creates the fibrin network that stops blood cells creating the clot. Finally fibroblasts migrate into the clot following chemotactic stimulation.

Therapeutic indications:

Chondral lesions treatment.

Contraindications:

Hypersensitivity to collagen.

Activity:

MeRG® carries out a "tent effect action" over mesenchymal cells, preventing their dispersion in the joint cavity. MeRG® is made of collagen fibres in a structure that favours cells adhesion. The three-dimensional structure of MeRG® enhances histological repair. In-vivo tests have shown that during the repair process, fibroblasts attach to the collagen fibrils, proliferate and orientate in order to reshape the damaged tissue. Collagen therefore supports tissue repair.

Degradation time:

MeRG® is physiologically degraded in 60/90 days. The fragments of collagen obtained by degradation are heat-sensitive and at a temperature of 37°C undergo a denaturation process, transforming into gelatine.

Product preparation:

Rehydrate the product with some drops of sterile physiological solution or with autologous bone marrow or PRP after shaping it.

MERG: COMPOSITION



The science of bone tissue

LABORATORIO R&D E ANALISI

10020 RIVA PRESSO CHIARI (TO) - Via G. Agnelli, 3 - Tel. +39 011 94 68 661 - Fax +39 011 94 64 036

ANALYSIS REPORT N.16219 - TA - 21/09/2016

RAPPORTO DI ANALISI N.16219 - TA - 21/09/2016

TITLE: *BIOCOLLAGEN COMPOSITION*

TITOLO: **COMPOSIZIONE BIOCOLLAGEN**

REFERENCE: *ASTM F2212 - 11*

RIFERIMENTO: ASTM F2212 - 11

INTERNAL PROCEDURE: *ILL - 007, ILL - 008, ILL - 014, ILL - 023*

RIFERIMENTO PROCEDURE INTERNE: ILL - 007, ILL - 008, ILL - 014, ILL - 023

TEST METHOD VALIDATION: *NO*

CONVALIDA METODO DI ANALISI: NO

REFERENCE DATA

DATI DI RIFERIMENTO

SAMPLES DESCRIPTION: *BIOCOLLAGEN MEMBRANE (10 mg)*

DESCRIZIONE CAMPIONI: MEMBRANA BIOCOLLAGEN (10 mg)

SAMPLES NUMBER: *3*

N° CAMPIONI: 3

PRODUCTION BATCH: *14290132*

LOTTO DI PRODUZIONE: 14290132

ANALYSIS PARAMETERS

PARAMETRI DI ANALISI

A) Proteins electrophoretic profile: No other bands than Collagen (As Type I Collagen Standard). B) GAGs spectrophotometric quantification ($\mu\text{g}/100\mu\text{g}$).
C) % Moisture.

ASSAYS	RESULTS
SAGGI ESEGUITI	ESITI
SENSITIVITY	
SENSIBILITA'	A) 0.01 ng/100mg B) 0.01 $\mu\text{g}/100\mu\text{g}$ C) 0.1 %
SAMPLES	
CAMPIONI	BIOCOLLAGEN 14290132 A) No other bands than Collagen are present.* B) 0.51 $\mu\text{g}/100\mu\text{g}$ * C) 7.2 % (%M)


*Reference A and B instruments output attached.

HEAD OF BIOCHEMISTRY LABORATORY
Responsabile del Laboratorio Biochimico
Dott. Paolo Fattori

The analysis report concern only the tested samples; the document can be partially reproduced only upon written approval by Bioteck Spa
Il Rapporto riguarda esclusivamente i campioni sottoposti a prova e non può essere riprodotto parzialmente salvo approvazione scritta da Bioteck Spa

MERG: GAGs QUANTIFICATION

Metodo: GAGs_QUANT_140916.mqa (520 nm)
 Modificati il: 14/09/2016 10:52:53 AM daWSSPECTRA\VERITON (SPECTRA)
 Spettrofotometro: GENESYS 6
 Numero di serie: 2M8F197001
 Insieme di..centrale :1.120
 Auto Zero: 14/09/2016 10:51:12 AM
 Misurato: 14/09/2016 10:54:14 AM daWSSPECTRA\VERITON (SPECTRA)
 Nome del f..risultati:BIOCOLLAGEN_GAGs.rqa



Standard

Numero	Concentrazione [µg/mL]	Ordinata [A]	Errore [A]	Già
1	0.00	0.006	0.006	Si
2	0.00	0.002	0.002	Si
3	1.00	0.024	0.002	Si
4	1.00	0.024	0.002	Si
5	5.00	0.112	0.001	Si
6	5.00	0.109	-0.002	Si

Campione	Diluizione Fattore	Ordinata [A]	Concentrazione [µg/mL]
SAMPLE 1	1	0.342	15.42
SAMPLE 2	20	0.057	51.41

MERG: PHYSIOCHEMICAL CHARACTERIZATION

Characteristic	Specification	Method	Unit
Chemical analysis			
Ashes (500-600° C)	< 5.0		%w
C	40.0-50.0	GC	%w
H	5.0-7.0	GC	%w
N	15.0-17.0	GC	%w
P	<1.0	GC	mg/g
Ca	<1.0	ICP-OES	mg/g
Mg	<0.5	ICP-OES	mg/g
Fe	<0.15	ICP	mg/g
Zn	<0.10	ICP	mg/g
Mn	<0.10	ICP	mg/g
Cu	<0.10	ICP	mg/g
Pb	<0.10	ICP	mg/g
As	<0.10	ICP	mg/g
Ca/P	Not appl.	Calculation	%
Cations and Anions			
Calcium (Ca ²⁺)	0.40-0.60	IC	mg/g
Magnesium (Mg ²⁺)	0.05-0.10	IC	mg/g
Sodium (Na ⁺)	0.30-0.70	IC	mg/g
Potassium (K ⁺)	<0.10	IC	mg/g
Ammonium (NH ₄ ⁺)	<0.20	IC	mg/g
Chloride (Cl ⁻)	<1.50	IC	mg/g
Nitrate (NO ₃ ⁻)	<0.10	IC	mg/g
Sulfate (SO ₄ ²⁻)	<0.30	IC	mg/g
Phosphate (PO ₄ ³⁻)	<0.05	IC	mg/g
Solution pH	5.0-6.0	IC	mg/g
Fat content			
	<0.01	GC/MS	%w
Collagen Characterization			
Purity of soluble collagen	No difference between bands of membrane derived collagen and standard	SDS-PAGE	
Elastin Assay	No elastin present	Western Blot	
Trypsin susceptibility		HPLC	
Amino acid analysis	Hydroxyproline >10.0 Tyrosine = 0.5 % Tryptophan = 0.5 %	HPLC	%w
Denatured collagen	=0.50%	HPLC	%w

Characteristic

Chemical analysis
Cations and Anions
Fat content
Collagen Characterization

Laboratory

Chemical Sciences Department
Instrumental Analysis Laboratory
University of Padua – Padua - Italy
Chelab S.r.l. – Resana (Treviso) – Italy

MERG CLINICAL WORK




Use of human mesenchymal stem cells (MSCs) combined to collagen-based scaffolds for cartilage and meniscus regeneration

Bighetti G, Sorbilli L, Faccini R., U.O. Ortopedia e Traumatologia Ospedale del Delta, Lagosanto, Ferrara, Italy

6 patients affected by intra-articular defects in the knee underwent cartilage reconstruction by implantation of MeRG® combined to autologous Bone Marrow Concentrate collected from the iliac crest.

At 6 and 12 months post-op, patients were evaluated from a clinical point of view, by using several scores (VAS, IKDC, Tegner/Lysholm); moreover, MRI analysis was also performed.

Bighetti G, Sorbilli L, Faccini R. Utilizzo di cellule mesenchimali autologhe associate a scaffold per la rigenerazione cartilaginea e meniscale. *OrthoAcademy, Tabloid Ortopedia, 2010; 17-18.*



← MeRG membrane combined to BMC and autologous thrombin

Clinical evaluation

At 12 months follow-up, any symptom affected the knee joint was nearly absent (VAS, IKDC and Tegner/Lysholm), patients began sports activity again, at levels comparable to pre-lesion condition.

MRI evaluation

From 6 to 12 months follow-up, defects were progressively filled with cartilaginous tissue showing structure comparable to healthy surrounding tissues; between cartilage and sub-chondral bone tissues a tide-mark like line was clearly visible.

Osteochondritis dissecans (OCD) of femoral condyle →



← 12 month post-op MRI evaluation of cartilage defect

Bighetti G, Sorbilli L, Faccini R. Utilizzo di cellule mesenchimali autologhe associate a scaffold per la rigenerazione cartilaginea e meniscale. *OrthoAcademy, Tabloid Ortopedia, 2010; 17-18.*


CLINICAL CASE REPORT



Arthroscopic Knee Cartilage Repair With Covered Microfracture and Bone Marrow Concentrate

Gigante A, Cecconi S, Ramazzotti D, Busilacchi A, Enea D. Polytechnic University of Marche, Orthopedics, Ancona, Italy
Calcagno S, Department of Orthopaedics, Sestri Levante Hospital (S. Calcagno), Sestri Levante, Italy.

37-yrs old man with medial joint-line pain in the left knee, MRI showed 3 cm² cartilage lesion on the medial femoral condyle. Cartilage defect was debrided, microfractures holes created, MeRG was combined to BM Concentrate, implanted and fixed in place with a mixture of fibrin glue and BMC. At 12 months post-op an MRI analysis was performed.




Clinical evaluation

At 12 months follow-up, any symptom affected the knee joint was absent. At 24 months, the patient was still asymptomatic.

MRI evaluation

MRI scan at 12 months showed good defect filling with a tissue signal very similar to that of surrounding tissue, as well as no signs of bone marrow edema.



Gigante A, Cecconi S, Calcagno S, Ramazzotti D, Busilacchi A, Enea D. Arthroscopic knee cartilage repair with covered microfracture and bone marrow concentrate. *Arthroscopy Techniques*, Vol 1, No 2 (December), 2012: pp e175-e180.

CLINICAL STUDY



Membrane-guided regeneration (MeRG) augmented with bone marrow concentrate (BMC) for cartilage repair in the knee. Histological results.

Enea D.¹, Calcagno S.², Alberto B.¹, Cecconi S.¹, Manzotti S.¹, Gigante A.¹

¹Polytechnic University of Marche, Orthopedics, Ancona, Italy, ²Rapallo Hospital, Orthopedics, Rapallo, Italy

5 consecutive patients affected by a focal isolated cartilage lesion in the knee underwent arthroscopic microfractures and the implant of MeRG[®] augmented with autologous BMC obtained from the iliac crest.



Arthroscopic evaluation

4 implants appeared nearly normal


1 implant appeared abnormal according the ICRS Cartilage Repair Assessment (CRA).

Results at the hystological evaluation

1 biopsy presented hyaline matrix (picture I)

1 biopsy presented a mixture of hyaline/fibrocartilage (picture J)

3 biopsies showed fibrocartilage (pictures F, G, H)



Each column represents a single patient.

Gigante A, Calcagno S, Cecconi S, Ramazzotti D, Manzotti S, Enea D. Use of collagen scaffold and autologous bone marrow concentrate as a one-step cartilage repair in the knee: Histological results of second-look biopsies at 1 year follow-up. *Int J Immunopathol Pharmacol* 2011;24:69-72.

CLINICAL STUDY




One-step cartilage repair in the knee: Collagen-covered microfracture and autologous bone marrow concentrate. A pilot study.

D. Enea ^{a,✉}, S. Cecconi ^a, S. Calcagno ^b, A. Busilacchi ^a, S. Manzotti ^a, A. Gigante ^a

(a) Department of Orthopedics, Polytechnic University of Marche, Via Tronto 10/A, 60020 Ancona, Italy

(b) Sestri Levante Hospital, Sestri Levante, GE, Italy

Nine patients with focal lesions of the condylar articular cartilage were consecutively treated with arthroscopic microfractures (MFX) covered with MeRG[®] membrane immersed in autologous bone marrow concentrate (BMC) from the iliac crest.




2nd look arthroscopy @ 1 year post-implant

Patients were retrospectively assessed using several standardized outcome assessment tools and MRI scans. Four patients consented to undergo second look arthroscopy and biopsy harvest.

Results at the histological evaluation

Cartilage macroscopic assessment at 12 months revealed that all the repairs appeared almost normal. Histological analysis showed a hyaline-like cartilage repair in one lesion, a fibrocartilaginous repair in two lesions and a mixture of both in one lesion



Biopsies stained with Safranin-O. Each column represents a single biopsy; lines e–h represent the chondral matrix; and lines i–l represent the osteochondral junction.


D. Enea, S. Cecconi, S. Calcagno, A. Busilacchi, S. Manzotti, A. Gigante. One-step cartilage repair in the knee: Collagen-covered microfracture and autologous bone marrow concentrate. A pilot study. *The Knee*, 2015;22:30–35.

CLINICAL EXPERIENCES (*P. Tessari, MD- Verona*)

2008-2013

45 patients (65% male, 35% female, 43 years mean age)

Knee	<i>medial compartment</i>	28	}	94%
	<i>lateral compartment</i>	5		
	<i>patella</i>			
Tibio-tarsal (talar area)		3	}	6%



Knee Medial C. Knee Lateral C. Patella Ankle

MeRG
Implantation sites

FOLLOW-UP AND RESULTS

Follow-up, 6 months post-op (45 patients - 100%)


<i>Excellent</i>		35	}	78%	Satisfying 90%
<i>Good</i>		5		12%	
<i>Fairly good</i>	4	8%	}	2%	Unsatisfying 10%
<i>Bad</i>		1			

Follow-up, 2 years post-op (42 patients - 94%)

<i>Excellent</i>		26	}	58%	Satisfying 83%
<i>Good</i>		11		25%	
<i>Fairly good</i>	5	11%	}	6%	Unsatisfying 13%
<i>Bad</i>		3			

CLINICAL IMAGES OF DEFECTS

A.W. - Male, 55 yrs




Medial femoral condyle


Defect with micro-fractures

A.C. - Male, 28 yrs
Football player

Iatrogenic calcinosis
of ligament



S.R. - Male, 57 yrs



B.R. - Female, 30 yrs, Tennis player



CLINICAL CASE - SECOND LOOK



Chondral lesion

Debridement of ill cartilage



Microfractures
months

MeRG positioning in CO₂

Second look at 7

CONSIDERATIONS

Dr. Tessari reported just 1 bad result (hyperemic reactive synovitis), assessed and arthroscopically treated .

At 2 years follow-up, he reported 3 bad results, just because 3 patients did not come back to controls.

However, results are aligned to what is reported in literature.

J. Gille, P. Behrens, P. Volpi, L. de Girolamo, E. Reiss, W. Zoch, and S. Anders. Outcome of Autologous Matrix Induced Chondrogenesis (AMIC) in cartilage knee surgery: data of the AMIC Registry. *Arch Orthop Trauma Surg.* 2013 January; 133(1): 87–93.

CONCLUSIONS

This 6 years-long experience with MeRG membrane has showed satisfying clinical outcomes in 83% of cases after 2 years.

All these cases had a complete pain resolution (VAS 0-2) and nearly complete joint recovery (ROM >95%).

Unsatisfying results were about 13% (6% just because referred to non controlled patients). All these patients had a story of soreness in the treated area and a subjective functional recovery that did not meet their expectations.

There was no case of infection. 4 synovitis, pharmacologically treated and in 2 cases arthroscopic revision and synovectomy were performed.

AMIC technique with MeRG membrane represents a *good and reproducible* surgical technique, in terms of timing and method.

The learning curve is quite fast before to easily manage the surgical procedure.

It is very important to strictly follow the indications, avoiding to extend the method to bigger defects or to cases associated with alteration of the loading axis and joint stability (ligament injuries).

Critical issues are micro-fractures of the sub-chondral bone which result in a very low supply of mesenchymal cells.

For this reason, we are thinking **to improve the quality of the regenerated tissue performing specific micro-perforations.**

NEW FINDINGS AND FUTURE TRENDS FOR MERG-BASED GCR (1/2)

Drilling and Microfracture Lead to Different Bone Structure and Necrosis during Bone-Marrow Stimulation for Cartilage Repair

Hongmei Chen,¹ Jun Sun,² Caroline D. Hoemann,¹ Viorica Lascau-Coman,¹ Wei Ouyang,¹ Marc D. McKee,³ Matthew S. Shive,² Michael D. Buschmann¹

¹Department of Chemical Engineering and Institute of Biomedical Engineering, Ecole Polytechnique de Montreal, P.O. 6079 Station Centre-ville, Montreal, QC, Canada H3C 3A7, ²BioSyntech Canada, Inc., Laval, Canada, ³Faculty of Dentistry, McGill University, Montreal, Canada

Drilling and Microfracture Lead to Different Bone Structure and Necrosis during Bone-Marrow Stimulation for Cartilage Repair

Hongmei Chen,¹ Jun Sun,² Caroline D. Hoemann,¹ Viorica Lascau-Coman,¹ Wei Ouyang,¹ Marc D. McKee,³ Matthew S. Shive,² Michael D. Buschmann¹

¹Department of Chemical Engineering and Institute of Biomedical Engineering, Ecole Polytechnique de Montreal, P.O. 6079 Station Centre-ville, Montreal, QC, Canada H3C 3A7, ²BioSyntech Canada, Inc., Laval, Canada, ³Faculty of Dentistry, McGill University, Montreal, Canada

Received 14 October 2009; accepted 23 March 2010
Published online 20 April 2010 in Wiley InterScience (www.interscience.wiley.com). DOI: 10.1002/jbm.b.32000

ABSTRACT: Bone marrow stimulation is performed using several surgical techniques that have and been anatomically compared or optimized for a desired cartilage repair outcome. In this study, we investigated acute osteoblast characteristics following microfracture and compared to drilling in a murine rabbit model of cartilage repair. Microfracture holes were made at a depth of 2 mm and drill holes to either 2 mm or 6 mm under cooled irrigation. Animals were sacrificed 1 day postoperatively and subchondral bone assessed by histology and micro-CT. The subchondral bone structure and bone density were compared between microfracture, microdrilling and drilling. Drilling holes were made and potentially sealed-off from the hole to provide access channels to marrow vessels. Our results showed that drilling could cause greater osteocyte death than microfracture due to increased pressure and osteoblasts, because osteocytes were present in contact with subchondral bone during drilling due to shearing and resulting of adjacent bone. Drilling deeper to 6 mm versus 2 mm prevented this effect and led to greater osteoblast formation. Our study revealed that differences between microfracture and drilling for an articular bone structure and osteocyte necrosis. Addition of cooling status suggest these differences significantly affect long-term cartilage repair outcomes. © 2010 Orthopaedic Research Society. Published by Wiley Periodicals, Inc. *J Orthop Res* 27: 1432–1438, 2009

Keywords: bone marrow stimulation; microfracture; Pringle drilling; subchondral bone; cartilage repair

A recent retrospective study on 25,128 knee arthroscopies demonstrated evidence of chondral defects in 13,074 knees (50%), and up to 67% were classified as localized and focal defects.¹ These cartilage defects have poor intrinsic healing capacity and, if left untreated, can lead to joint degeneration, chronic pain, and disability. Several surgical repair methods have been developed to treat focal cartilage defects.^{2–7} The most widely practiced are bone marrow stimulation techniques such as subchondral drilling,⁸ autologous arthroplasty,⁹ and microfracture (MF).¹⁰ The rationale behind these methods is to initiate bleeding and wound repair via small fractures in subchondral bone. Drilling, first introduced by Pridie in 1959,⁸ was prevalent for decades¹¹ prior to MF becoming advocated in the 1990s.^{12–14} Both of these procedures aim at carrying bone marrow elements to repair cartilage defects.^{15,16} MF uses an arthroscopic tool to penetrate the subchondral plate, and is thought to be superior to Pringle drilling,^{10,17,18} which employs a hand-driven or motorized drill that may produce heat necrosis in the bone—a hypothesis that has never been tested. The MF technique represents a first-line treatment option for full-thickness cartilage defects in the knee,¹⁹ showing a decrease in pain and improved knee function after 1- or 2-year follow-up.^{20–22} Recent clinical data also provided evidence that outcomes after MF are similar to those of autologous chondrocyte transplantation.^{23–25}



by the formation of a blood clot. The subsequent classical wound repair cascade comprised of an acute inflammatory response and cell chemotaxis leads to the generation of a vascularized granulation tissue, and the proliferation of chondrogenic mesenchymal progenitor cells with a capacity to differentiate into multiple morphological cell types.²⁶ These remodeling can proceed along with the initiation of chondrogenesis that resembles endochondral ossification, to form new bone in deeper zones and fibrous or fibrocartilaginous tissues in the more superficial chondral regions.^{10,15,18} A number of studies using bone marrow stimulation have been carried out to improve cartilage repair by modifying the perforation technique,²⁷ combining with periosteal grafts,²⁸ hybrid implants,^{29–32} scaffolds,³³ or other surgical methods.³⁴ However, to our knowledge, no previous study has specifically investigated the characteristics of the acute fracture themselves, such as hole shape, extent of hematoma, extent of bone damage, morphology of fractured bone and bone marrow, and osteocyte necrosis. Furthermore, no controlled animal or clinical study has directly compared MF to drilling techniques despite their clinical popularity.

The aim of the current study was to compare acute osteoblast characteristics 24 h after microfracture (MF) or microdrilling (MD) procedures in a laboratory murine rabbit model. We tested the hypothesis that 1) microfracture induces bone compaction around holes, and 2) microdrilling induces bone necrosis around holes. The effect of hole depth was also examined where MF was performed to a depth of 2 mm (MF2), while MD was performed to a depth of either 2 mm (MD2) or 6 mm (MD6). Our results revealed significant differences between MF and MD, which could potentially influence subsequent repair response and long-term cartilage repair properties.

Correspondence to: Michael D. Buschmann (E-mail: michael.buschmann@epm.polymtl.ca)
J Orthop Res 27: 1432–1438, 2009. © 2010 Orthopaedic Research Society. Published by Wiley Periodicals, Inc.

1432

This study reveals “significant differences between Microfracture (MF) and Microdrilling (MD), which could potentially influence subsequent repair responses and longer term cartilage repair properties”.




“Microfracture with an awl induced fracturing and bone compaction around holes that were largely sealed-off from adjacent bone marrow, in contrast to drilling which cleanly removed bone debris and left channels that communicate between the hole and marrow. Microfracture also produced a high level of osteocyte necrosis in adjacent bone, in contrast to our drilling method which included cooled irrigation and did not cause apparent thermal necrosis.”

MD 6mm >> MD 2mm > MF 2mm

MERG REHABILITATION PROTOCOL


- Immobilization with brace and unloading for 10 days
- Passive motion with articulated brace and unloading for 10 days
- Progressive loading with use of crutches for 10 days, associated to active joint mobilization
- Complete loading and joint exercises to achieve full ROM, for 3 weeks
- Restarting of sports activity, 3 months after surgery




- Immobilization with fixed brace and unloading for 1 week
- Mobilization with articulated brace and progressive load with crutches for 2 weeks
- Complete loading with muscle strengthening exercises, achieving full ROM in 1 week
- Return to normal activity after 1 month (no sports activity)
- Return to sports activity, after 2 months



This protocol is indicated for non load-bearing areas, as well as for regions of indirect loading (e.g. patella)



MERG: CRITICAL ISSUES IN REHAB




PERICARDIUM MEMBRANES

Heart Patch


Acellular native collagen patch from equine pericardium.

Soft, strong and elastic due to its structure of the fibers oriented in all directions.




SEM pictures show the compact multi-layer structure and the dense network of collagen fibers.

SEM pictures by Dept. Of Biology, University of Padova, Italy




The use of Heart Patch in the reinforcement of the posterior tibialis tendon.


HEART® PATCH INDICATION




HEART PATCH IN POSTERIOR TIBIALIS TENDON REPAIR




The patch is cut to the needed size, before rehydration



The pericardium strip is passed under the tendon



The pericardium strip is sutured and the excess is then trimmed



The tendon is completely wrapped by the pericardium strip

HEAD AND NECK GRAFT LINE

BIOOTECK
The science of bone tissue

Osteoplast Spine
the biological choice

Enzymatic deantigenation
preserves bone biological properties of bone

A unique, innovative deantigenation process

Heart DM
natural re-absorbable membrane

Application of Heart DM
in repairing a dural defect

Biochemical characterisation

Osteoplast Spine
A complete line of grafts
for Spine surgery

BIOOTECK
The science of bone tissue

Heart DM
biological matrix for dural replacement

A natural membrane with high biomechanical performance


Application of Heart DM in repairing a dural defect

Biochemical characterisation

Heart DM
Equine Pericardium
membrane for dura
substitution and repair




Heart DM
Available in 5 different size from
30x20mm to 120 mmx x160



Heart DM as dura substitute
in a meningioma case














HEART DM INDICATIONS



CODE AND AVAILABLE SIZES

<u>Code</u>	<u>Description</u>
HRT-40DM	HEART® PERICARDIUM MEMBRANE 25 X 30 mm
HRT-41DM	HEART® PERICARDIUM MEMBRANE 50 X 50 mm
HRT-42DM	HEART® PERICARDIUM MEMBRANE 60 X 80 mm
HRT-43DM	HEART® PERICARDIUM MEMBRANE 80 X 140 mm
HRT-44DM	HEART® PERICARDIUM MEMBRANE 120 X 160 mm


TABLE OF MAIN INDICATIONS

	SMALL BONE LOSS FRACTURES	NON-UNIONS MAL-UNIONS	ACETABULAR REVISIONS HIP REVISIONS	SPINE FUSIONS JOINT FUSIONS	FEMORAL TIBIAL OSTEOTOMY	ARTICULAR CARTILAGE DEFECTS	ROTATOR CUFF TENDONS REINFORCEMENT DURA REPAIR
	●	●	●				
	●						
			●				
			●				
	●		●				
				●	●		
		●		●			
	●			●			
				●			
				●			
	1 ● 2	1 ● 2	1 ●	● 2			
							● ●
						●	

CLINICAL CASES 1

BONE SUBSTITUTES IN ACETABULAR REVISION

Courtesy of prof. Roberto Sessa, MD. – Catania, Italy



Acetabular cup revision in a 74 years old female – RX pre-op.



Osteoplast cancellous chips .
RX 3 months post-op.




Osteoplast cancellous chips .
RX 18 months post-op.


CLINICAL CASES 2

BONE SUBSTITUTES IN ACETABULAR REVISION

Courtesy of prof. Roberto Sessa, MD. – Catania, Italy




OTC-C8 Flex Cortical Sheet.




OSP-070 Flex Acetabular Mat

In presence of acetabular protusio, flex cortical sheet is used as a new acetabular wall associated to the acetabular mat for filling the defect



RX 2 months post-op.



RX 8 months post-op.




RX 16 months post-op.

G. Sessa et al.: *Equine bone tissue in acetabular revision: our experience*
MINERVA ORTOP TRAUMATOL 2010;61:469-76)

CLINICAL SERIE

BONE SUBSTITUTES IN ACETABULAR REVISION

Courtesy of Prof. Massimiliano Marcucci
Director of C.E.S.A.T. Centro Eccellenza Sostituzioni Articolari Fucecchio (FI) - Italy



RX pre-op

***Bone regeneration in revision hip arthroplasty using equine-derived bone grafts: a retrospective study**

Massimiliano Marcucci¹, Angelo Graceffa¹, Nicola Piolanti¹, Pier Francesco Indelli² and Leonardo Latella¹


¹Articular Replacements Excellence Center (CESAT) - Fondazione Onlus "...In Cammino...", Fucecchio, Italy

²Stanford University, San Jose, California, USA


**Paper submitted for publication*

Abstract


Background: During the last decade, total hip arthroplasty has become a common procedure performed in younger patients, as well as elderly ones. This has led to an increase in total hip arthroplasty revisions. Loosening of primary components with associated bone loss represents the major cause of total hip arthroplasty revision. This study evaluates the safety and performance of an enzyme-deantigenic equine-derived bone graft material when used as an alternative to bone autografts and/or allografts in acetabular defect reconstruction. **Methods:** Records of 55 patients who were treated for Paprosky Type II or III acetabular bone defects with arthroplasty revisions using equine-derived bone and followed for an average of 18 months (range 6–48 months) were analyzed. **Results:** Of the 59 revisions, 53 (89.8%) were regarded as successful, with evidence of incorporation of the equine bone graft. Failures included one infection (1.7% of revisions) and five cases of aseptic loosening of the acetabular prosthesis (8.5% of revisions). These results are consistent with those of studies having a similar follow-up period for allografts used in combination with trabecular metal components. **Conclusion:** Results of the present study suggest that enzyme-treated equine-derived bone grafts may be a valid alternative to autogenous and homologous bone grafts in total hip arthroplasty revision.



Acetabular prosthesis to revise



Curettage of the defect



Reaming of the acetabulum




Acetabulum ready for grafting



*Activagen Injectable Paste
mixed with Osteoplast
chips*




Filling of the acetabular defect




The defect is filled




*Chips are compacted by using
the reamer with opposite
rotation*



*Appearance of the reconstruct
acetabular cavity*




Acetabular prosthesis is placed




The defect is filled



Chips are compacted by using the reamer with opposite rotation



Appearance of the reconstruct acetabular cavity



Acetabular prosthesis is placed



RX pre-op.



RX post-op.

CLINICAL CASES

BONE SUBSTITUTES IN HIP REVISION

Courtesy of prof. Roberto Sessa, MD. – Catania, Italy



RX pre-op.



Impaction grafting - RX post-op.



RX 30 months post-op.

CLINICAL CASE

ANEURISMATIC BONE CYST IN A 9 Y. O. CHILD

**Equine-derived bone substitute in orthopedics and traumatology:
authors' experience**


S. SANTINI. P. BARBERA. M, MODENA. R. SCHIAVON. M. BONATO



X-Ray pre-op.



Part of the graft is still visible
X-ray 6 months post-op.




Perfect integration of the graft.
RX 20 months post-op.


CLINICAL CASE

INFECTED SEVERE NON-UNION IN UPPER LIMB


Pictures are courtesy of Dr. F. Da Rin - Codivilla Putti Hospital, Cortina d'Ampezzo (BL)



Bioteck chips mixed with autografts and bone marrow



Antibiotic spacer is removed. The mixture have been put in a cut syringe for some minutes in order to get a more solid graft due to blood coagulation and then placed under the metal plate.



Graft positioning



X-ray control 6 months post-op.

CLINICAL SERIES

PROXIMAL HUMEROUS MULTI-FRAGMENTS FRACTURE

M. Trono, , G.Lucidi

U.O. Ortopedia e Traumatologia - Ospedale Infermi, Rimini

*Abstract

Background: This retrospective study aimed to evaluate the functional outcomes and incidence of complications in cases of proximal humeral fractures treated by open reduction and internal fixation with locking plates and concomitant equine-derived bone graft.

Methods: One hundred seventeen patients showing Neer 2-4 proximal humeral fractures were treated with locking plate fixation and concomitant grafting with an equine-derived bone substitute. Clinical and radiologic assessments were performed at the 1, 3, and 12-month follow-up appointments to evaluate fracture healing, alignment and reduction, and occurrence of complications. Functional recovery was assessed at the 12-month follow-up with the Constant-Murley scoring system.

Results: Radiologic and clinical bone healing was achieved in all patients. The overall complication rate was 5.1% and included nonunions, postoperative bleeding, and mobilization of the osteosynthesis devices. At the 12-month follow-up, the mean Constant score was 87.7 ± 8.5 and the mean Functional Recovery, in comparison with the contralateral unaffected arm, was $94.8\% \pm 4.2$.

Conclusions: The use of equine bone blocks may be highly advisable in conjunction with locking plate fixation for the treatment of proximal humeral multipart fractures.

Level of Evidence: This work is a case series investigating the results of a treatment and corresponding to a Level IV of Evidence.



73 years old female - RX pre-op



RX 6 months post-op



6 months clinical result

* Paper submitted for publications


CLINICAL CASE 1

SEVERE NON-UNION IN LOWER LIMBS

Courtesy of dr. Bruno Di Maggio, MD. – Piedimonte Matese (CE), Italy




Bifocal femur fracture in a 32 years old woman of 120 kgs.



Osteoplast chips + PRP

Previous treatment in an other hospital
had no result after one year



Implanted materials:
Osteoimplant chips + autologousPRP




RX post-op



RX 4 months post-op



Clinical control 8 months post-op




RX 15 months post-op



CLINICAL CASE 2

SEVERE NON-UNION IN LOWER LIMBS


Courtesy of dr. Bruno Di Maggio, MD. – Piedimonte Matese (CE), Italy



Double open leg fracture in a 20 years old woman.
External fixation in urgency, but signs of bone infection after 3 month



Implanted materials:
Osteoimplant chips + autologousPRP



RX post-op



RX 3 months post-op



Clinical control 6 months post-op



CLINICAL CASE

HTO WITH EQUINE BONE WEDGE

Courtesy of Rainero Del Din, MD. – SS. Antonio e Margherita Hospital, Tortona



Pre-op RX



The wedge is adjust to final
shape and size
Bone marrow concentrated
cells are used to enrich it



The wedge is placed into the
Tibial osteotomic lacuna



A metal plate is positioned and fixed with screws



Post-op Xray



2 months post-op Xray



4 months post-op Xray



CLINICAL SERIES

*A novel equine-derived pericardium membrane for dural repair: a preliminary, short-term investigation.

Roberto Centonze MD, Emiliano Agostini MD, Samantha Massaccesi MD, Stefano Toninelli MD, Letterio Morabito MD

Division of Neurosurgery, Azienda Ospedaliera Marche Nord, Pesaro, Italy


Background: A large variety of biological and artificial materials are employed in dural repair, each of them with major limitations: autologous grafts have limited availability and involves an additional incision and surgical time; cadaveric preparations and heterologous materials entail the risk of *iatrogenic transmission of prions* whereas synthetic substitutes have been reported to cause inflammatory reactions and graft rejection. An equine-derived pericardium membrane has been developed (Heart^o Bioteck, Vicenza, Italy) with mechanical and safety-related features that would make it suitable for neurosurgical application.

Aims: This preliminary study aimed to evaluate the short-term safety and efficacy of Heart^o membrane in dural repair procedures following meningioma surgeries.

Methods: The medical records of 8 patients surgically treated for an intracranial meningioma and undergoing duraplasty with the Heart membrane were reviewed retrospectively. Clinical and radiological assessments were performed on the 1st and 30th post-operative days. The occurrence of graft-related complications, such as cerebrospinal fluid (CSF) leakage, post-operative hematoma, wound infections, meningitis and neurological symptoms were analysed.

Results: A watertight closure was achieved in all the patients. Post-operatively, no patients exhibited CSF leak, cerebral contusion, haemorrhage or wound infections. The one-month follow-up revealed no evidence of pseudomeningocele, wound breakdown or meningitis. *Neurologic complications were observed in 3 patients but not directly imputable to the dural substitute or its application*

Conclusions: The pericardium membrane allowed to achieve watertight dural closure without graft-related adverse events in all the patients. Further investigations should be performed to assess medium- and long-term clinical outcomes in a larger set of patients.



TC pre-op



The membrane is sutured for covering the defect




TC post-op


* Paper submitted for publications

CLINICAL CASE HEART PATCH 1 AVULSION OF TRICEPS BRACHII


Courtesy of Andrea Atzei, MD. – Hand Surgeon, Treviso





Bone Anchor fixation



Heart® Patch to bridge the gap




Heart® Patch secured to Bone Anchor



Wrapping the Triceps Rafe Tendon





Core Sutures cross Heart® Patch



Soft tissue wrap around the Patch

18 months P.O.

Full ROM
90% strength



CLINICAL CASE HEART PATCH 2 AUGMENTATION OF CHRONIC RUPTURE DISTAL BICEPS TENDON

Courtesy of Andrea Atzei, MD. – Hand Surgeon, Treviso


male - 28 yo heavy
laborer

Closed injury palmar
aspect right dominant
elbow




misdiagnosed as
muscular "contusion" of
the biceps brachii

left untreated.

After 2 months
severe impairment of
elbow flexion




Sub-total rupture of biceps brachii at teno-muscular junction 5 cm gap




Two core sutures (Ethibond #2)



Preparation of the graft cutting Heart® Patch



Augmentation with Heart® Patch



Ethibond 2-0 interrupted sutures at periphery and central



16 months P.O.

Full ROM
90% strength

CLINICAL CASE HEART PATCH 3 AUGMENTATION OF ANULAR LIGAMENT RECONSTRUCTION


Courtesy of Andrea Atzei, MD. – Hand Surgeon, Treviso

20 y-o female

13 years earlier
post-traumatic isolated
Proximal Radio-Ulnar
Dislocation





Anular Ligament







Lateral Approach
To PRU-j




Radial Head
Resection & Reduction

Ligament remnants








Anterior Approach to PRU-j





Harvesting bone tunnel
for Graft Fixation





Positioning Graft around Radial Neck







Suturing Graft To ligament remnants with interrupted sutures

Bone Anchor to fix Graft To the bone tunnel





Suturing ligament
with Graft
And securing
with anchors sutures



Pre-op.



6 months p.o.



CLINICAL CASE HEART PATCH 4 INTERPOSITION ARTHROPLASTY POST-TRAUMATIC DISPLASTIC RADIAL HEAD RESECTION

Courtesy of Andrea Atzei, MD. – Hand Surgeon, Treviso


male - 16 yo
aerospatial student
(aircraft pilot)

Radial head fracture
left non-dominant elbow


early arthritis
of the
capitellum





Restricted elbow motion



Radial head Resection





Resurfacing /interposition Arthroplasty

18 months P.O.

ROM 90%

90% strength

Stable PRUJ

